

EBP06485 - IMMUNIZING PEPTIDE

SKU: EBP06485



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_003684.2; NP_663325.1; NP_663327.2
Immunogen	Peptide with sequence SELDPKGQHVC, from the internal region of the protein sequence according to NP_003684.2; NP_663325.1; NP_663327.2.
Product Comments	The peptide was used in the production of EB06485 - Goat Anti-SCARF1 / SREC Antibody. It is currently untested in blocking applications.
Peptide Sequence	SELDPKGQHVC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)