

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06493 - IMMUNIZING PEPTIDE

SKU: EBP06493

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

ID

NP_055416.2

 $\textbf{Immunogen} \begin{tabular}{l} \textbf{Peptide with sequence CRLVPPSKRRHKGSA, from the C Terminus of the protein sequence according to NP_055416.2.} \end{tabular}$

Product The peptide was used in the production of EB06493 - Goat Anti-EHD2 Antibody. It is currently untested in blocking

Comments applications.

Peptide

Sequence

CRLVPPSKRRHKGSA

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade