

EBP06507-B - IMMUNIZING PEPTIDE

SKU: EBP06507-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001440.2; NP_001153171.1
Immunogen	Peptide with sequence CNKRFFVHNEQVY., from the internal region of the protein sequence according to NP_001440.2; NP_001153171.1.
Product	The peptide was used in the production of EB06507-B - Goat Anti-FHL1 / SLIM1, Biotinylated Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	CNKRFVHNEQVY.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)