

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06522 - IMMUNIZING PEPTIDE

SKU: EBP06522

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession

 $NP_005684.2; \ NP_009052.4; NP_001123573.1; \ NP_001238863.1; \ NP_001350524.1$ ID

Peptide with sequence CRLQDKKLPPLLSEI, from the internal region of the protein sequence according to

Immunogen NP_005684.2; NP_009052.4;NP_001123573.1; NP_001238863.1; NP_001350524.1.

The peptide was used in the production of EB06522 - Goat Anti-LXR alpha / LXR beta Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

CRLQDKKLPPLLSEI Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade