

# EBP06532 - IMMUNIZING PEPTIDE

SKU: EBP06532



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_853549.1; NP\_874387.1; NP\_874388.1; NP\_874393.1; NP\_874394.1; NP\_877570.1; NP\_877571.1; NP\_877572.1; NP\_877573.1; NP\_878270.1; NP\_898883.1

**Immunogen** Peptide with sequence LETLKDICSPKTDY, from the C Terminus of the protein sequence according to NP\_853549.1; NP\_874387.1; NP\_874388.1; NP\_874393.1; NP\_874394.1; NP\_877570.1; NP\_877571.1; NP\_877572.1; NP\_877573.1; NP\_878270.1; NP\_898883.1.

**Product Comments** The peptide was used in the production of EB06532 - Goat Anti-CREM Antibody. It is currently untested in blocking applications.

**Peptide Sequence** LETLKDICSPKTDY

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

## DOCUMENTS

- [Data Sheet](#)