

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06532 - IMMUNIZING PEPTIDE

SKU: EBP06532

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

NP_853549.1; NP_874387.1; NP_874388.1; NP_874393.1; NP_874394.1; NP_877570.1; NP_877571.1; Accession NP_877572.1; NP_877573.1; NP_878270.1; NP_898883.1 ID

Peptide with sequence LETLKDICSPKTDY, from the C Terminus of the protein sequence according to

Immunogen NP 853549.1; NP 874387.1; NP 874388.1; NP 874393.1; NP 874394.1; NP 877570.1; NP 877571.1;

NP 877572.1; NP 877573.1; NP 878270.1; NP 898883.1.

The peptide was used in the production of EB06532 - Goat Anti-CREM Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

LETLKDICSPKTDY Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade