

EBP06539 - IMMUNIZING PEPTIDE

SKU: EBP06539



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_060121.3; NP_001128304.1; NP_001337433.1
Immunogen	Peptide with sequence PTAESEAKVKTC, from the N Terminus of the protein sequence according to NP_060121.3; NP_001128304.1; NP_001337433.1.
Product	The peptide was used in the production of EB06539 - Goat Anti-AHI1 / Jouberein Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	PTAESEAKVKTC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)