

EBP06541 - IMMUNIZING PEPTIDE

SKU: EBP06541



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_060621

Immunogen Peptide with sequence CNSESDSSEAKEEG, from the internal region of the protein sequence according to NP_060621.

Product Comments The peptide was used in the production of EB06541 - Goat Anti-RIF1 (aa 1551-1563) Antibody. It is currently untested in blocking applications.

Peptide Sequence CNSESDSSEAKEEG

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png