

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06560 - IMMUNIZING PEPTIDE

SKU: EBP06560

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_057051; NP_036543.4; NP_001180360.1; NP_001295226.1 ID

 $\label{eq:local_$

The peptide was used in the production of EB06560 - Goat Anti-LYP / PTPN22 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

CPPNKPAESVQSNNS Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade