

EBP06563 - IMMUNIZING PEPTIDE

SKU: EBP06563



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_054772.1
Immunogen	Peptide with sequence C-KTVMLSKQSESAI, from the C Terminus of the protein sequence according to NP_054772.1.
Product	The peptide was used in the production of EB06563 - Goat Anti-FLVCR Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-KTVMLSKQSESAI
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)