

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06571 - IMMUNIZING PEPTIDE

SKU: EBP06571

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_031394.2 ID

Immunogen

Peptide with sequence C-QYKRDKFKKTKYGSQ, from the internal region (near the C Terminus) of the protein

sequence according to NP_031394.2.

The peptide was used in the production of EB06571 - Goat Anti-RASA3 / GAP1IP4BP Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

C-QYKRDKFKKTKYGSQ Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade