

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06572 - IMMUNIZING PEPTIDE

SKU: EBP06572

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_002799.3; NP_001265637.1; NP_001265638.1 ID

Peptide with sequence C-VILRKNPNYDL, from the C Terminus of the protein sequence according to NP_002799.3; Immunogen NP_001265637.1; NP_001265638.1.

The peptide was used in the production of EB06572 - Goat Anti-TRAP2 / Proteasome subunit 26S Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-VILRKNPNYDL Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade