

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06595 - IMMUNIZING PEPTIDE

SKU: EBP06595

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 000745.1; NP 001128633.1; NP 001128634.1; NP 009294.1 ID

Peptide with sequence GDTKEQRILNHVLQC, from the N Terminus of the protein sequence according to Immunogen

NP_000745.1; NP_001128633.1; NP_001128634.1; NP_009294.1.

The peptide was used in the production of EB06595 - Goat Anti-COMT (N Terminus) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

GDTKEQRILNHVLQC Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade