

EBP06595 - IMMUNIZING PEPTIDE

SKU: EBP06595



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000745.1; NP_001128633.1; NP_001128634.1; NP_009294.1
Immunogen	Peptide with sequence GDTKEQRILNHVLQC, from the N Terminus of the protein sequence according to NP_000745.1; NP_001128633.1; NP_001128634.1; NP_009294.1.
Product	The peptide was used in the production of EBP06595 - Goat Anti-COMT (N Terminus) Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	GDTKEQRILNHVLQC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)