

EBP06636 - IMMUNIZING PEPTIDE

SKU: EBP06636



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_066264.3; NP_065574.3; NP_001136405.1
Immunogen	Peptide with sequence CKEKATRPSQGHQP, from the C Terminus of the protein sequence according to NP_066264.3; NP_065574.3; NP_001136405.1.
Product	The peptide was used in the production of EBP06636 - Goat Anti-Choline acetyltransferase Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	CKEKATRPSQGHQP
Shipping	
Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)