

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06643 - IMMUNIZING PEPTIDE**

**SKU:** EBP06643

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_000356.1; NP\_001152759.1; NP\_001244955.1 ID

 $Peptide \ with \ sequence \ C-LKPEFVDIINAKQ, \ from \ the \ C \ Terminus \ of \ the \ protein \ sequence \ according \ to \ NP\_000356.1;$ Immunogen NP\_001152759.1; NP\_001244955.1.

The peptide was used in the production of EB06643 - Goat Anti-Triosephosphate isomerase Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

C-LKPEFVDIINAKQ Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite\_plus\_medium.png$ Grade