

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06645 - IMMUNIZING PEPTIDE

SKU: EBP06645

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_001074.2; NP_246273.2; NP_236914.2; (NP_237223) ID

Peptide with sequence C-LINGESGQAKRN, from the C Terminus of the protein sequence according to NP_001074.2; Immunogen NP_246273.2; NP_236914.2; (NP_237223).

The peptide was used in the production of EB06645 - Goat Anti-PDE5A Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

C-LINGESGQAKRN Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade