

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06659 - IMMUNIZING PEPTIDE

SKU: EBP06659

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_000656.1; NP_001269378.1 ID

Peptide with sequence QFDHYSKQDRCSDL, from the C Terminus of the protein sequence according to NP_000656.1;

The peptide was used in the production of EB06659 - Goat Anti-ACHE Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

QFDHYSKQDRCSDL Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade