

EBP06659 - IMMUNIZING PEPTIDE

SKU: EBP06659



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000656.1; NP_001269378.1
Immunogen	Peptide with sequence QFDHYSKQDRCSL, from the C Terminus of the protein sequence according to NP_000656.1; NP_001269378.1.
Product Comments	The peptide was used in the production of EB06659 - Goat Anti-ACHE Antibody. It is currently untested in blocking applications.
Peptide Sequence	QFDHYSKQDRCSL
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)