

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06675 - IMMUNIZING PEPTIDE**

**SKU:** EBP06675

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

ID

NP\_000368.1

Immunogen Peptide with sequence C-SPADKKRSGKKKI, from the internal region of the protein sequence according to NP\_000368.1.

**Product** The peptide was used in the production of EB06675 - Goat Anti-Wiskott-Aldrich Syndrome / WASP Antibody. It is

**Comments** currently untested in blocking applications.

**Peptide** 

Sequence

C-SPADKKRSGKKKI

**Shipping** 

Ambient Instructions

**Product** 

Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite\_plus\_medium.png$