

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06684 - IMMUNIZING PEPTIDE**

**SKU:** EBP06684

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Shipped at ambient temperature, store at -20°C.

**Accession**  $NP\_001341935.1; \ NP\_954654.1; \ NP\_001341936.1; \ NP\_001341939.1$ ID

Peptide with sequence C-QEAIQDLWQWRKSL, from the C Terminus of the protein sequence according to

Immunogen NP\_001341935.1; NP\_954654.1; NP\_001341936.1; NP\_001341939.1.

The peptide was used in the production of EB06684 - Goat Anti-Nucleophosmin / NPM1 Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

C-QEAIQDLWQWRKSL Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade