

EBP06684 - IMMUNIZING PEPTIDE

SKU: EBP06684



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001341935.1; NP_954654.1; NP_001341936.1; NP_001341939.1
Immunogen	Peptide with sequence C-QEAIQDLWQWRKSL, from the C Terminus of the protein sequence according to NP_001341935.1; NP_954654.1; NP_001341936.1; NP_001341939.1.
Product Comments	The peptide was used in the production of EB06684 - Goat Anti-Nucleophosmin / NPM1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-QEAIQDLWQWRKSL
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)