

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06695 - IMMUNIZING PEPTIDE

SKU: EBP06695

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession

 $NP_003992.3;\ NP_001002273.1;\ NP_001002274.1;\ NP_001002275.1;\ NP_001177757.1$ ID

Peptide with sequence C-PDALEEPDDQNRI, from the C Terminus of the protein sequence according to NP_003992.3; Immunogen NP_001002273.1; NP_001002274.1; NP_001002275.1; NP_001177757.1.

The peptide was used in the production of EB06695 - Goat Anti-CD32 / FCGR2B Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-PDALEEPDDQNRI Sequence

Shipping

Instructions Ambient

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade