

EBP06696 - IMMUNIZING PEPTIDE

SKU: EBP06696



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_733821.1; NP_733822.1
Immunogen	Peptide with sequence PRTQSPQNCSIM, from the C Terminus of the protein sequence according to NP_733821.1; NP_733822.1.
Product	The peptide was used in the production of EB06696 - Goat Anti-Lamin A / Lamin A/C Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	PRTQSPQNCSIM
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png