

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06696 - IMMUNIZING PEPTIDE

SKU: EBP06696

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 733821.1; NP 733822.1 ID

Peptide with sequence PRTQSPQNCSIM, from the C Terminus of the protein sequence according to

The peptide was used in the production of EB06696 - Goat Anti-Lamin A / Lamin A/C Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide PRTQSPQNCSIM

Sequence Shipping

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade