

EBP06725 - IMMUNIZING PEPTIDE

SKU: EBP06725



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001129.1; NP_015531.1
Immunogen	Peptide with sequence CRKLG TAMNPCSRT, from the C Terminus of the protein sequence according to NP_001129.1; NP_015531.1.
Product Comments	The peptide was used in the production of EB06725 - Goat Anti-AGRP / Agouti-related protein Antibody. It is currently untested in blocking applications.
Peptide Sequence	CRKLG TAMNPCSRT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)