

EBP06745 - IMMUNIZING PEPTIDE

SKU: EBP06745



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004304.1; NP_945355.1; NP_001244257.1; NP_001244259.1; NP_001244260.1; NP_001316993.1
Immunogen	Peptide with sequence C-HDHIPLRPQS, from the internal region of the protein sequence according to NP_004304.1; NP_945355.1; NP_001244257.1; NP_001244259.1; NP_001244260.1; NP_001316993.1.
Product Comments	The peptide was used in the production of EB06745 - Goat Anti-Arrestin beta 2 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-HDHIPLRPQS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)