

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06751 - IMMUNIZING PEPTIDE

SKU: EBP06751

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

ID

NP_002737.2; NP_001035145.1; NP_001103361.1

Immunogen Peptide with sequence GGEPRRTEGVGP-C, from the N Terminus of the protein sequence according to NP_002737.2; NP_001035145.1; NP_001103361.1.

The peptide was used in the production of EB06751 - Goat Anti-MAPK3/ERK1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

Sequence

GGEPRRTEGVGP-C

Shipping

Ambient Instructions

Product

Grade

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$