

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06755 - IMMUNIZING PEPTIDE**

**SKU:** EBP06755

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

100 μg **Unit Size** 

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 079496.1; NP 001269376.1 ID

Peptide with sequence CQSRRCRKNAFQEL, from the internal region of the protein sequence according to Immunogen

NP 079496.1; NP 001269376.1.

The peptide was used in the production of EB06755 - Goat Anti-ADAM33 Antibody. It is currently untested in **Product** 

Comments blocking applications.

**Peptide** 

CQSRRCRKNAFQEL Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade