

EBP06803 - IMMUNIZING PEPTIDE

SKU: EBP06803



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000367.1; NP_001017535.1; NP_001351014.1;NP_001017536.1
Immunogen	Peptide with sequence CGNQDYKYRVSD, from the internal region of the protein sequence according to NP_000367.1; NP_001017535.1; NP_001351014.1;NP_001017536.1.
Product	The peptide was used in the production of EB06803 - Goat Anti-VDR Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	CGNQDYKYRVSD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)