

EBP06816 - IMMUNIZING PEPTIDE

SKU: EBP06816



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001899; NP_680090; NP_680091; NP_680092; NP_680093
Immunogen	Peptide with sequence C-DELVNYVNRNT, from the internal region of the protein sequence according to NP_001899; NP_680090; NP_680091; NP_680092; NP_680093.
Product Comments	The peptide was used in the production of EB06816 - Goat Anti-Cathepsin B Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-DELVNYVNRNT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png