

EBP06823 - IMMUNIZING PEPTIDE

SKU: EBP06823



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_003168.2; NP_001128524.1
Immunogen	Peptide with sequence C-KYLQQNRHVKDKN, from the internal region of the protein sequence according to NP_003168.2; NP_001128524.1.
Product	The peptide was used in the production of EB06823 - Goat Anti-Spleen tyrosine kinase / SYK Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-KYLQQNRHVKDKN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)