

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06831 - IMMUNIZING PEPTIDE**

**SKU:** EBP06831

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

100 μg **Unit Size** 

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 055211.2; NP 542198.2 ID

Peptide with sequence C-DHLKKDSQNKTLLK, from the internal region of the protein sequence according to Immunogen

NP\_055211.2; NP\_542198.2.

The peptide was used in the production of EB06831 - Goat Anti-VPS41 Antibody. It is currently untested in **Product** 

**Comments** blocking applications. **Peptide** 

Sequence

C-DHLKKDSQNKTLLK

Shipping

Grade

Instructions Ambient

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png