

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06859 - IMMUNIZING PEPTIDE

SKU: EBP06859

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP_665906.1; NP_665908.1; NP_002122.1; NP_665909.1; NP_665912.1; NP_665910.1 ID

Peptide with sequence C-LASKQEKDGTEK, from the internal region of the protein sequence according to Immunogen

NP 665906.1; NP 665908.1; NP 002122.1; NP 665909.1; NP 665912.1; NP 665910.1.

The peptide was used in the production of EB06859 - Goat Anti-HMGA1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-LASKQEKDGTEK Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png