

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06860 - IMMUNIZING PEPTIDE

SKU: EBP06860

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

ID

NP_079141.2

Peptide with sequence C-QNQRMKSKRWQKNN, from the internal region of the protein sequence according to Immunogen NP_079141.2.

The peptide was used in the production of EB06860 - Goat Anti-NANOG Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

Sequence

C-QNQRMKSKRWQKNN

Shipping

Ambient Instructions

Product

Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$