

EBP06876 - IMMUNIZING PEPTIDE

SKU: EBP06876



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_057446.1; NP_001128413.1; NP_001128416.1; NP_001128417.1; NP_001128418.1
Immunogen	Peptide with sequence C-QDILWEEAKEAPADK, from the C Terminus of the protein sequence according to NP_057446.1; NP_001128413.1; NP_001128416.1; NP_001128417.1; NP_001128418.1.
Product Comments	The peptide was used in the production of EBP06876 - Goat Anti-Ghrelin preproprotein Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-QDILWEEAKEAPADK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)