

EBP06877 - IMMUNIZING PEPTIDE

SKU: EBP06877



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000234.1; NP_001185465.1
Immunogen	Peptide with sequence C-EHLKCLRKSGEEQ, from the internal region of the protein sequence according to NP_000234.1; NP_001185465.1.
Product Comments	The peptide was used in the production of EBP06877 - Goat Anti-MEFV Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-EHLKCLRKSGEEQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)