

EBP06878 - IMMUNIZING PEPTIDE

SKU: EBP06878



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001104491.1; NP_001104492.1; NP_001104493.1; NP_033951.2
Immunogen	Peptide with sequence C-PSQSDVVETE, from the internal region of the protein sequence according to NP_001104491.1; NP_001104492.1; NP_001104493.1; NP_033951.2.
Product Comments	The peptide was used in the production of EB06878 - Goat Anti-Runx1 (mouse) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PSQSDVVETE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)