

EBP06904 - IMMUNIZING PEPTIDE

SKU: EBP06904



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Unit Size | 100 µg |
| Storage Instructions | Shipped at ambient temperature, store at -20°C. |
| Accession ID | NP_000852.1; NP_001091681.1; NP_001091682.1; NP_001091683.1 |
| Immunogen | Peptide with sequence CNENFKKTFKRILH, from the C Terminus of the protein sequence according to NP_000852.1; NP_001091681.1; NP_001091682.1; NP_001091683.1. |
| Product Comments | The peptide was used in the production of EB06904 - Goat Anti-Histamine Receptor H1 (C Term) Antibody. It is currently untested in blocking applications. |
| Peptide Sequence | CNENFKKTFKRILH |
| Shipping Instructions | Ambient |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png |

DOCUMENTS

- [Data Sheet](#)