

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06904 - IMMUNIZING PEPTIDE**

**SKU:** EBP06904

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Shipped at ambient temperature, store at -20°C.

**Accession** 

 $NP\_000852.1; \ NP\_001091681.1; \ NP\_001091682.1; \ NP\_001091683.1$ ID

Peptide with sequence CNENFKKTFKRILH, from the C Terminus of the protein sequence according to NP\_000852.1; Immunogen NP\_001091681.1; NP\_001091682.1; NP\_001091683.1.

The peptide was used in the production of EB06904 - Goat Anti-Histamine Receptor H1 (C Term) Antibody. It is **Product** 

**Comments** currently untested in blocking applications.

**Peptide** 

**CNENFKKTFKRILH** Sequence

Shipping Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade