

EBP06940-B - IMMUNIZING PEPTIDE

SKU: EBP06940-B



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_001116423.1; NP_002027.2

Immunogen Peptide with sequence HRAELSPSTENSSQLDFEDC., from the N Terminus of the protein sequence according to NP_001116423.1; NP_002027.2.

Product Comments The peptide was used in the production of EB06940-B - Goat Anti-Duffy / FY / DARC, Biotinylated Antibody. It is currently untested in blocking applications.

Peptide Sequence HRAELSPSTENSSQLDFEDC.

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)