

EBP06940 - IMMUNIZING PEPTIDE

SKU: EBP06940



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001116423.1; NP_002027.2
Immunogen	Peptide with sequence HRAELSPSTENSSQLDFED-C, from the N Terminus of the protein sequence according to NP_001116423.1; NP_002027.2.
Product Comments	The peptide was used in the production of EB06940 - Goat Anti-Duffy / FY / DARC Antibody. It is currently untested in blocking applications.
Peptide Sequence	HRAELSPSTENSSQLDFED-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)