

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP06944 - IMMUNIZING PEPTIDE**

**SKU:** EBP06944

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_004554.3; NP\_001278485.1 ID

 $\label{eq:lmmunogen} \textbf{Immunogen} \begin{tabular}{ll} \textbf{Peptide with sequence KPWKPGDKEPCAH, from the internal region of the protein sequence according to NP\_004554.3; NP\_001278485.1. \end{tabular}$ 

**Product** The peptide was used in the production of EB06944 - Goat Anti-PCK2 / PEPCK-M Antibody. It is currently untested in

**Comments** blocking applications.

**Peptide** 

**KPWKPGDKEPCAH** Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite\_plus\_medium.png$ Grade