

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06945 - IMMUNIZING PEPTIDE

SKU: EBP06945

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 057167.2; NP 001153732.1; NP 001153731.1; NP 001153730.1; NP 001153698.1; NP 149421.1 ID

Peptide with sequence C-HTSEDGKVQVTRPDQ, from the internal region of the protein sequence according to Immunogen

NP 057167.2; NP 001153732.1; NP 001153731.1; NP 001153730.1; NP 001153698.1; NP 149421.1.

The peptide was used in the production of EB06945 - Goat Anti-Cannabinoid Receptor 1 Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-HTSEDGKVQVTRPDQ Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png