

EBP06945 - IMMUNIZING PEPTIDE

SKU: EBP06945



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_057167.2; NP_001153732.1; NP_001153731.1; NP_001153730.1; NP_001153698.1; NP_149421.1
Immunogen	Peptide with sequence C-HTSEdGKVQVTRPDQ, from the internal region of the protein sequence according to NP_057167.2; NP_001153732.1; NP_001153731.1; NP_001153730.1; NP_001153698.1; NP_149421.1.
Product Comments	The peptide was used in the production of EB06945 - Goat Anti-Cannabinoid Receptor 1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-HTSEdGKVQVTRPDQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)