

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP06975 - IMMUNIZING PEPTIDE

SKU: EBP06975

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_065162.1; NP_851994.1; NP_001186776.1 ID

Immunogen Peptide with sequence C-PSEGSLEDEELPAQ, from the internal region of the protein sequence according to NP_065162.1; NP_851994.1; NP_001186776.1.

The peptide was used in the production of EB06975 - Goat Anti-SMURF1 Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

C-PSEGSLEDEELPAQ Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade