

EBP06986 - IMMUNIZING PEPTIDE

SKU: EBP06986



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_149107.4
Immunogen	Peptide with sequence C-KEVDLTRLSEKERRD, from the internal region of the protein sequence according to NP_149107.4.
Product Comments	The peptide was used in the production of EB06986 - Goat Anti-NEK9 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KEVDLTRLSEKERRD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)