

EBP06993 - IMMUNIZING PEPTIDE

SKU: EBP06993



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_004091.3; NP_742103.1; NP_742101.2; NP_001287941.1; NP_001287942.1
Immunogen	Peptide with sequence EDSQDLNEQSVKKTC, from the N Terminus of the protein sequence according to NP_004091.3; NP_742103.1; NP_742101.2; NP_001287941.1; NP_001287942.1.
Product	The peptide was used in the production of EB06993 - Goat Anti-EYA4 Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	EDSQDLNEQSVKKTC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)