

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07025 - IMMUNIZING PEPTIDE

SKU: EBP07025

SPECIFICATIONS

Formulation $100\mu g$ of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

NP_000255.1; NP_001077072.1; NP_001077071.1; NP_001077073.1; NP_001077074.1; NP_001077075.1; Accession

NP_001077076.1 ID

Peptide with sequence ELQDVECEERPR, from the C Terminus of the protein sequence according to

Immunogen NP 000255.1; NP 001077072.1; NP 001077071.1; NP 001077073.1; NP 001077074.1; NP 001077075.1;

The peptide was used in the production of EB07025 - Goat Anti-PTCH (C Terminus) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

ELQDVECEERPR Sequence

Shipping

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade