

EBP07040 - IMMUNIZING PEPTIDE

SKU: EBP07040



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_073753.3
Immunogen	Peptide with sequence C-PRQPKKFQKGRKDN, from the internal region of the protein sequence according to NP_073753.3.
Product Comments	The peptide was used in the production of EBP07040 - Goat Anti-CDCP1 (isoform 1, internal) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PRQPKKFQKGRKDN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png